

Human Cryopreserved Hepatocytes

Lot number: HUM4064A

Date: December 30, 2013



Cryo Characterization Report (CCR)


Lot Overview		
Qualification	Catalog Number	Quantity
Cryopreserved human hepatocytes, suspension qualified	HUCSD	5.0x10 ⁶ viable cells/vial

Donor Demographics								
Sex	Race	Age	BMI	Tobacco Use	Alcohol Use	Drug Use	Serological Data	Cause of Death
Female	Caucasian	48	36	Yes	Yes	No	Negative	Stroke

Additional donor demographic information, including relevant medical and medication history, is available upon request

Post-thaw Viability and Cell Quality Assessment			
Thawing Medium Used	Centrifuge Conditions	% Viability (post-thaw)	Viable cell yield per vial
Cryopreserved Human Hepatocyte Thawing Medium	100 x g for 8 min	80%	5.0x10 ⁶

HepatoMeter report available upon request



HUM4064A, Post-Thaw, 10X

Human Cryopreserved Hepatocytes

Lot number: HUM4064A

Date: December 30, 2013



<i>Suspension Metabolism</i>			
Isoform	Substrate	Marker Metabolite	pmol/min/10 ⁶ cells
<i>General</i>	<i>100 μM 7-Ethoxycoumarin</i>	<i>7-Hydroxycoumarin</i>	<i>838.0</i>
<i>SULT</i>	<i>100 μM 7-Hydroxycoumarin</i>	<i>7-Hydroxycoumarin Sulfate</i>	<i>29.2</i>
<i>UGT</i>	<i>100 μM 7-Hydroxycoumarin</i>	<i>7-Hydroxycoumarin Glucuronide</i>	<i>1981</i>

Contact customer service to place an order or to obtain additional information on any of our lots. This may include supplementary donor demographic information, current inventory, and photomicrographs at multiple timepoints and magnifications.

To contact TRL:

Main: 800-748-8979 or 919-549-3580

Sales: 919-549-3593

Email: customerservice@trlcells.com

Web: www.triangleresearchlabs.com