

Human Cryopreserved Hepatocytes
 Lot number: HUM4145
 Date: 12/17/15



Cryo Characterization Report (CCR)

Lot Overview	
Qualification	Catalog Number
Cryopreserved human hepatocytes, Induction	HUCPI

Storage Conditions: <-150°C


Donor Demographics								
Sex	Race	Age	BMI	Tobacco Use	Alcohol Use	Drug Use	Serological Data	Cause of Death
Male	Caucasian	1 yrs 8 mos	16.4	No	No	No	Negative	Anoxia 2nd Drowning

Additional donor demographic information, including relevant medical and medication history, is available upon request


Post-thaw Viability and Cell Quality Assessment			
Thawing Medium Used	Optimal Centrifuge Conditions	% Viability (post-thaw)	Viable cell yield per vial
MCHT50	100 g x 8 minutes	98%	5.8 x 10 ⁶

Monolayer Assessment					
Plating Medium Used	Well Format	Culture Medium Used	Optimal Seeding Density	Initial Attachment Efficiency	Monolayer Confluency @ 120 hrs
MP250	24 well	MM250	0.8x10 ⁶ /mL	90%	100%

Characterization was completed in a 24-well collagen-coated plate with extra-cellular matrix overlay



HUM4145, 72hrs, 20X



HUM4145, 96hrs, 10X

<i>Induction</i>			
Isoforms	Control Inducer	Fold Induction Specific Activity	Fold Induction mRNA Expression
<i>CYP1A2</i>	<i>50µM Omeprazole</i>	<i>29.2</i>	<i>337.0</i>
<i>CYP2B6</i>	<i>1mM Phenobarbital</i>	<i>14.1</i>	<i>18.8</i>
<i>CYP3A4</i>	<i>10µM Rifampicin</i>	<i>9.2</i>	<i>13.5</i>

Cryopreserved human hepatocytes were thawed and plated on 24-well collagen I coated plates, overlaid with Matrigel©, then dosed in triplicate with vehicle control (0.1% DMSO) or control inducers for 72 hours. The fold induction was calculated by dividing the induced level by the vehicle control level.

Media products used for characterization:

MCHT50 - Cryopreserved human hepatocyte thawing medium, 50mL

MM250 - Hepatocyte maintenance medium, 250mL

MP250 - Human hepatocyte plating medium, 250mL

Contact customer service to place an order or to obtain additional information on any of our lots. This may include supplementary donor demographic information, current inventory, and photomicrographs at multiple timepoints and magnifications.

To contact TRL:

Main: 800-748-8979 or 919-549-3580

Sales: 919-549-3593

Email: customerservice@trlcells.com

Web: www.triangleresearchlabs.com