

Human Cryopreserved Hepatocytes

Lot number: HUM4232

Date (D/M/Y): 20/06/2017



Cryo Characterization Report (CCR)


Lot Overview		
Qualification	Catalog Number	Manufacture Date (D/M/Y)
Cryopreserved human hepatocytes, General Plateable	HUCPG	09-04-2017

Storage Conditions: <-150°C


Donor Demographics								
Sex	Race	Age	BMI	Tobacco Use	Heavy Alcohol Use	Drug Use	Serological Data	Cause of Death
Male	Caucasian	32 yrs	21	No	No	Yes	Negative	Drug Overdose

Additional donor demographic information, including relevant medical and medication history, is available upon request

Post-thaw Viability and Cell Quality Assessment			
Thawing Medium Used	Optimal Centrifuge Conditions	% Viability (post-thaw)	Viable cell yield per vial
MCHT50	100 g x 8 minutes	82	10 x 10 ⁶



HUM4230A, 24hrs, 10X



HUM4230, 96hrs, 20X

Media products used for characterization:

MCHT50 - Cryopreserved human hepatocyte thawing medium, 50mL

MCHT50P - Cryopreserved pooled human hepatocyte thawing medium, 50mL

Contact customer service to place an order or to obtain additional information on any of our lots. This may include supplementary donor demographic information, current inventory, and photomicrographs at multiple timepoints and magnifications.

To contact Lonza:

Customer Service: +1 800-638-8174 | order.us@lonza.com

Scientific Support: +1 800-521-0390 | scientific.support@lonza.com | Web: www.lonza.com

EU/Row Scientific Support: +32 87 321 611 | scientific.support.eu@lonza.com;

Lonza.com/TRLheps