


## Cryo Characterization Report (CCR)

Lot Overview			
Qualification	Catalog Number	Passage	Quantity
Cryopreserved human liver Stellate cells	HUCLS	P0	1.0x10 <sup>6</sup> viable cells/vial

Donor Demographics								
Sex	Race	Age	BMI	Tobacco Use	Alcohol Use	Drug Use	Serological Data	Cause of Death
Male	Caucasian	26	22	Yes; unknown product/frequency	None	None	Negative	Cardiac Arrest

*Additional donor demographic information, including relevant medical and medication history, is available upon request*

*Each lot of primary cells is tested found non-reactive to HIV-1, HIV-2, HTLV I, HTLV II, Hepatitis B, Hepatitis C Syphilis and CMV. However, no known test can offer complete assurance that these infectious agents are not present. Since we cannot test all the pathogens, please treat the culture as a potential infectious reagent at Biosafety Level 2 or higher. Human hepatic non-parenchymal cells are cultured in DMEM/F12 containing 10% FBS and antibiotics. FACs analysis is shown below for P0 culture:*



Contact customer service to place an order or to obtain additional information on any of our lots. This may include supplementary donor demographic information, current inventory, and photomicrographs at multiple timepoints and magnifications.

**To contact TRL:**

Main: 800-748-8979 or 919-549-3580  
 Sales: 919-549-3593  
 Email: customerservice-trl@lonza.com  
 Web: www.trlcells.com