

Cryopreserved Rat (Wistar Han) Hepatocytes -  
Plateable

Catalog Number: RWCP01

Lot number: RWH139

Date: 07/11/16



**LONZA**

Triangle Research Labs – Now Part of Lonza

FOR RESEARCH USE ONLY. CAUTION: Not intended for diagnostic or therapeutic use. Treat all animal derived products as potentially pathogenic. Wear personal protective equipment and practice appropriate laboratory and disposal techniques for potentially pathogenic or bio-hazardous material.

<i>Donor Demographics</i>			
Species	Strain	Sex	Age
<i>Rat</i>	<i>Wistar Han</i>	<i>Male</i>	<i>5 weeks</i>


Storage Conditions: <-150°C

<i>Post-thaw Viability and Cell Quality Assessment</i>			
Thawing Medium Used	Optimal Centrifuge Conditions	% Viability	Viable cell yield per vial
<i>MCRT50</i>	<i>100g x 10min</i>	<i>87%</i>	<i>6.8 x 10<sup>6</sup></i>

<i>Monolayer Assessment</i>				
Plating Medium Used	Well Format	Culture Medium Used	Initial Attachment Efficiency	Monolayer Confluency @ 48 hrs
<i>Hepatocyte Plating Medium (MP100)</i>	<i>24 well</i>	<i>Hepatocyte Maintenance Medium (MM250)</i>	<i>80%</i>	<i>80%</i>

<i>Suspension Metabolism</i>			
Isoform	Substrate	Marker Metabolite/ *Substrate Clearance	Activity (pmol/min/10 <sup>6</sup> cells)
<i>General</i>	<i>100 uM 7-Ethoxycoumarin</i>	<i>7-Ethoxycoumarin Clearance*</i>	<i>N/A</i>
<i>General</i>	<i>100 μM 7-Ethoxycoumarin</i>	<i>7-Hydroxycoumarin</i>	<i>10.9</i>
<i>SULT</i>	<i>100 uM 7-Hydroxycoumarin</i>	<i>7-Hydroxycoumarin Sulfate</i>	<i>350.9</i>
<i>UGT</i>	<i>100 uM 7-Hydroxycoumarin</i>	<i>7-Hydroxycoumarin Glucuronide</i>	<i>721.6</i>

All incubations were concluded at 30 minutes






RWH139, 24hrs, 10X

RWH139, 48hrs, 10X

Triangle Research Labs ▲ 6 Davis Drive, Research Triangle Park, NC 27709 ▲

1.919.549.3593

To place an order or inquire about our products and services, contact us by phone:  
919.549.3593  or  
Email: [customerservice-trl@lonza.com](mailto:customerservice-trl@lonza.com); or visit us on the web at [www.trlcells.com](http://www.trlcells.com).

Triangle Research Labs  6 Davis Drive, Research Triangle Park, NC 27709   
1.919.549.3593 